

## Biocidal ZF™

**N-86620**

Use biocides carefully. Always read the label and product information before use.

### Spray Disinfectant for Incubators and Sterile Cabinets in Cell Culture Area

#### Cell Culture Area



#### Effectiveness:

Biocidal ZF™ is active against bacteria, fungi, spores and enveloped viruses.

#### Non-volatile:

Active ingredients are non-volatile. They protect cell cultures from microbial contamination and do not invade cell cultures via air. Thus cell cultures are protected against contamination and the disinfectant itself.

#### Application for incubators

- Empty incubator
- Spray Biocidal ZF™ on all surfaces and let settle
- Put cultures right back

Surfaces will be protected for 2-4 weeks. (Before disinfecting you may clean incubator with Biocidal ZF™ if necessary)

#### Application in sterile benches

- Empty sterile bench
- Then clean with Biocidal ZF™ and finally spray on all surfaces.

Protection will prevail until next workshift. No combustible gases can form, unlike by use of ethanol.

Classification according to EU-Regulation No. 1272/2008: Harmful to aquatic life with long lasting effects.

<b>Active ingredients:</b>	Quaternary ammonium compounds -free of aldehydes-		
<b>Storage:</b>	Room temperature		
<b>Materials compatibility:</b>	Stainless steel		
<b>Cat.#</b>	WAK-ZF-1	1 x 1	liter spray bottle
	WAK-ZF-6	6 x 1	liter spray bottle
	WAK-ZF-12	12 x 0.2	liter spray bottle

L-01-001-02-0-01

Siemensstraße 9  
D-61449 Steinbach/Ts. • Germany  
Tel. +49 61 71 / 28 43 0 • Fax +49 61 71 / 69 56 46  
e-mail: info@wak-chemie.com  
www.wak-chemie.com



## WAK - Chemie Medical GmbH

